



# EPIDURAL ANESTHESIA NEEDLES

## “FIXED WING”

Features :

- The curved design of the needle point allows protection of the catheter from damage during procedures and smooth insertion into the epidural space.
- The ultrathin wall of the cannula makes it possible to insert a large catheter.
- Gradual insertion at 1cm intervals allows confirmation of the insertion depth.
- Arrows on the fixed wing assist in awareness of needle placement.
- It's plastic stylet is fitted to prevent tissue coring and is lighter to enhance tactile feel during insertion.



GAUGE	LENGTH(mm)	COLOUR
16G	90mm	WHITE
17G	90mm	VIOLET
18G	90mm	PINK

 **Dr. Japan Co., Ltd.**

1-8-5 Shinjuku, Shinjuku-ku, Tokyo, Japan

<http://www.drjapan-jp.com/>