


▶ Dr. Japan Face Shield



We want to protect everyone from COVID-19
with **Japanese Highest level Technology !!**



It is used very often in Japanese medical institutions
Because...

1. Transparency

Manufactured from highly transparent A-PET material,
it is **98.82% transparent**

2. Anti-fog coating

The film is coated to prevent it from
becoming cloudy

3. Splash prevention

Wide range of shield designed based on
splash evaluation results



*Medical field needed
a more protective face shield*

Made in Japan

*Developed a face shield with Akamatsu Kasei,
a leading food packaging company in Japan.*

Product Specifications


Name	Dr. Japan Face Shield
Item No	AK-003 Dr.Japan-SP
Quantity	120 pc (24pc×5bags) / export carton
Material	A-PET (Amorphous polyethylene terephthalate) Heatproof temperature 60°C
Size	Shield 0.30mm×355mm×260mm Belt 35mm×748mm
Weight	38.7g/pcs
Carton	478mm×188mm×279mm 6.6Kg

Product Concept

- ◆ Shield size that wraps the entire face ensures protection
- ◆ Outstanding transparency reduces eye strain and improves work efficiency
- ◆ Anti-fog processing is applied, so you can always work with high transparency
- ◆ Easy to store and dispose of integrally molded PET material
- ◆ You can work with glasses and goggles

※ Only the face side (inner surface) is anti-fog processed. Please note that wiping with alcohol will reduce the anti-fog effect
 ※ It may be deformed by heat. Please avoid distribution, storage, and its use in an environment of 60 degrees or higher

How to Use



Manufacturer

AKAMATSU KASEIKOUGYOU Co.,Ltd.

119-1, Mitsuho-kaitaku, Mitsuho-aza, Matsushige-cho,
Itano-gun, Tokushima 771-0298
TEL +81-88-699-3733 FAX +81-88-699-3732

Distributor

Dr. Japan Co., Ltd.

1-1 Kagurazaka, Shinjuku-ku, Tokyo 162-0825
TEL +81-3-3513-8766 FAX +81-3--3513-8226
info.overseas@drjapan-jp.com